




LWave[®] – Möglichkeiten und Anwendungsbeispiele der laserakustischen Oberflächenwellen -Spektroskopie

—
Dr. Stefan Makowski
Fraunhofer Institute for Material and Beam Technology IWS, Germany



Contents

Introduction

Measurement Principle, Evaluation Strategies, the Role of Young's Modulus for Characterization

Applications in the Thin Film Technology

Nanometer Films

Structural Changes in Thin Films

Subsurface Damage of Silicon Wafers

Quality Control Concept for Superhard Carbon Coatings

Pores in Metal Films

Low-K Films

Nitriding Depth - Life Time of Ammonia Reactors

AI-Assisted Identification of Different Types of Coating Failures


Methodical Aspects, Development

Comparison with Nanoidentation; Current Development





Introduction



Laser-induced Surface Acoustic Wave Spectroscopy



Wikipedia: "Seismic waves", public domain





„Miniature earthquake machine“



- Young's Modulus
 - Coating thickness
 - Density
 - Depth of: Nitriding layer, Case-hardening, Damage layers
 - Porosity, crack density
 - Delamination
 - Build-up structure
- Anything that affects your mechanical integrity of the material

Laser-induced surface acoustic wave spectroscopy - LAwave




- (1) Short laser pulse
→ broadband surface acoustic waves
- (2) SAW propagation,
depth depends on frequency
- (3) SAW measurement: piezoelectric element
- (4) Fourier transformation
→ dispersion curve ($c = f(f)$) + other data



LAwave setup for small parts
© J. Jeibmann/ Fraunhofer IWS

Effective Young's modulus of heterogeneous materials



Defects reducing Young's modulus

- Porosity (thermal spraying, AM)
- Integrity, cracks (damage layers, bad adhesion)
- Structure (columnar growth, 3D printing)



Applications

Application - Overview

Coatings




- Thin coatings: PVD, CVD, ALD, spin coating, electroplating, ...
- Thick coatings: thermal-spraying, laser-cladding, ...
- Metal films
- Polymeric sensor films

Semiconductors

- Wafer bonding layers; damage depth
- Low-k films

Bulk and surfaces

- Case hardening, nitriding
- Rolling, hammering, shot peening
- 3D printing
- Composition (e.g. carbide content in hard metal)



Case study: Very thin films < 10 nm

Material

- PVD coatings with thickness < 10 nm

Results


- Measurement of Young's Modulus

▪ HfO ₂	220.4 GPa
▪ Native SiO ₂	39.8 GPa
▪ SiO _x	41.7 GPa
▪ a-C	373.4 GPa
▪ TiAlN	142.8 GPa
▪ Silicon wafer	165.2 GPa (C11)
▪ Silicon wafer (high doping)	162.9 GPa (C11)

- Measurement of thickness of Si/Al/Al₂O₃ multilayer stack


▪ Thermal Al ₂ O ₃	3.9 nm
--	--------

➔ Possibilities beyond conventional nanoindentation




Case study: Structural changes in thin films

Growth dependent density of RF-sputtered refractory metals



T.A. Yakupov, et al., Materials Today: Proceedings 4 (2017) 4469–4476.

Influence of element concentration on Young's modulus in AlYN wurtzite thin films



D. I. Solonenko, et al., Applied Physics A 131,690, (2025)



Case study: Subsurface damage in semiconductor wafers

Material

- Semi-conductor surfaces, damaged from processing

Results

- Damage layer \rightarrow dispersion
- Slope = damage layer depth \rightarrow allows quantification



Case study: Quality control of superhard carbon coatings


Material

- Superhard amorphous carbon coatings (ta-C, H-free DLC), hardness 40..70 GPa
- Application: Low-wear low-friction coating, e.g. piston pins in ICE, motorcycle chain
- State-of-the art: Nanoindentation → slow and error-prone technique with high indenter wear


Results

- LAwave allows to access
 - Coating modulus, coating hardness
 - Coating thickness

in less than 60 seconds



Case study: Quality control of superhard carbon coatings II



Lw ta-C on piston pin

- Position left
 - Analysis
 - Model
- Position center
 - Analysis
 - Model
- Position right
 - Analysis
 - Model

+Compare Dispersion curve - Comparison

MATERIAL MODEL

Films	ρ / g/cm ³	E / GPa	ν	t / μ m
Diamond like Carbon	3,00	545,12	0,19	2,12

Substrate	ρ / g/cm ³	E / GPa	ν
PFluegel	7,61	204,20	0,29


Bounds

6,33 MHz

83,97 MHz

- QUALITY CONTROL**
- Coating Thickness
 - Coating Young's Modulus


Material Model



→ EXPORT MODEL TO CLIPBOARD

Example: Pores in metal films

Material: 2 μm Ti coating on Si wafer, deposited with plasma activated electron beam deposition



➔ **Effective Young's Modulus strongly correlates with porosity observed in SEM cross section**


Example: Porous low-k films

Material

- Nano-porous SiCOH low-k films
- High porosity: > 40 %
- Rel. permittivity $k < 2.5$
- Minimum required stiffness $E > 5$ GPa

Results

- Young's modulus and density can be measured
- Higher reliability than results from nanoindentation



Irradiation with 172 nm and 222 nm photons, $\rho = 1.2 \text{ gcm}^{-3}$ and $d = 200 \text{ nm}$

Prager et al. *Microelectronic Engineering* 85 (2008) 2094–2097



Case study: Nitriding depth in NH₃ reactor

Material

- Steel wall of ammonia reactor
- Nitriding depth depends on processing time and conditions


Application

- Nitriding depth measurement
- Study nitride hardening depth
- Life-time, e.g. in ammonia exposed reactors



Outlook: Identification of defects

WC-Co by thermal spraying





Methodical Aspects


Comparison with instrumented indentation testing

Coating Materials


- Porous low-k films

Result


- Effective modulus is strongly overestimated with indentation due to compressed pores



Surface acoustic waves
Reversible deformation
→ True elasticity



Indentation
Densification of microdefects
→ Distorted results



Comparison with instrumented indentation testing (nanoindentation)

	LWave	Nanoindentation
Method	Dynamic: Sound velocity $c \sim \sqrt{E/\rho}$	Quasi-static: $E_r \sim dP/dh$
Measuring area	$> 5 \times 5 \text{ mm}^2$ (integral method)	$< 10 \text{ }\mu\text{m}^2$ (local method)
Measuring time	$t \sim 1 \text{ min}$	$t \gg 1 \text{ min}$ (incl. preparation and calibration)
Minimal film thickness	$d > 1 \text{ nm}$	$d > 100 \text{ nanometers}$
Surface roughness	No requirements	Smooth surface necessary
Difficult material systems	Transparent and high damping materials; coated foils	Soft and superhard materials, very thin coatings
Benefits	Thin and porous coatings	True plastic behavior; high lateral resolution

Development and Technology Readiness

State of the art



- Standalone system, manual handling, quality control concepts - TRL 9
- 40+ systems world wide

New features under development


- Automated measurement – TRL4
- In-situ measurement up to 600 °C – TRL 4
- AI assisted evaluation – TRL 5
- Mobile head for robot or hand measurement

Fraunhofer IWS provides

- 35+ years of experience in LAwave application across all industries
- Development partnership for industry and academia
- Job measurements, feasibility studies, quality control concepts
- Custom solutions for research, quality control, analysis and automation



Summary



LWave®

- Mechanical properties of coatings and surfaces
- Fast and non-destructive
- Coatings from 1 nm to 1 mm
- Used for R&D and Quality Control
- 35+ years of experience; 40+ systems world-wide, 70+ peer reviewed contributions



LAwave around the world

Distribution USA

Fraunhofer USA
44792 Helm Street
Plymouth, MI 48170 USA
Phone: +1 734-354-4335
lawave@fraunhofer.org

Headquarter Germany

Fraunhofer Institute for Material and
Beam Technology IWS
Winterbergstrasse 28, 01277 Dresden,
Germany
Dr. Martin Zawischa
+49 351 83391 3096
lawave@iws.fraunhofer.de

Distribution Japan

Sanyo Trading Co., Ltd.
2-11, Kanda Nishiki-cho, Chiyoda-ku Tokyo
101-0054 Japan
Koya Tamori
+81-3-3518-1196
k-tamori@sanyo-trading.co.jp

Distribution China

GermanTech Nano Co. Ltd.
Room 1408, Tower B, Techart Plaza
No.30 Xue Yuan Road, Haidian District,
100083 Beijing, P.R.China
+86 10 82867920/21/22
contact@germantech.com.cn

Image: Aris Katsaris CC BY-SA 3.0

Contact

lawave@iws.fraunhofer.de

Dr. Martin Zawischa +49 351 83391 3096
Dr. Stefan Makowski +49 351 83391 3192

Fraunhofer Institute for Material and Beam Technology IWS
Winterbergstraße 28
01277 Dresden
Germany

Links and Resources

LAWave Website: [Link](#)

One page product sheet: [Link](#)

Two minute video on Youtube: [Link](#)

Scientific Literature using LAWave: [Link](#)